

Secure your care

Samsung Healthcare Cybersecurity

Bringing peace of mind to your hospital and patients

To address the emerging need for cybersecurity, Samsung provides a solution to support our customers by offering the tools to protect against cyberthreats that may compromise invaluable patient data and ultimately degrade the quality of care. Samsung's Cybersecurity Solution strives to abide by the CIA triad (Confidentiality, Integrity, and Availability) and takes a comprehensive approach to providing impeccable protection with the following pillars: Intrusion prevention, Access control, and Data protection.



Intrusion prevention

- Tools for protecting against cyber threats from external attacks
- Security tools include Anti-virus & Firewall
 - Secured operating system



Access control

- Strengthened surveillance for tracking the access of patient information
- Account management
 - Enhanced audit trail



Data protection

- Encryption functions for safeguarding data whether at-rest or in-transit
- Data protection
 - Transmission security

About Samsung Medison CO., LTD.

Samsung Medison, an affiliate of Samsung Electronics, is a global medical company founded in 1985. With a mission to bring health and well-being to people's lives, the company manufactures diagnostic ultrasound systems around the world across various medical fields. Samsung Medison has commercialized the Live 3D technology in 2001 and since being part of Samsung Electronics in 2011, it is integrating IT, image processing, semiconductor and communication technologies into ultrasound devices for efficient and confident diagnosis.

- This product, features, options and transducers are not commercially available in all countries.
- Sales and shipments are effective only after the approval by the regulatory affairs. Please contact your local sales representative for further details.
- This product is a medical device, please read the user manual carefully before use.
- S-Vue Transducer™ is the name of Samsung's advanced transducer technology.
- Optical Disk Drive is not available for this product.

1. Optional feature which may require additional purchase.
2. Clinical image acquired by the HS60 V1.00 ultrasound system.
3. Clinical image acquired by the HS60 V2.00 ultrasound system.
4. Clinical image acquired by the HS60 V2.02 ultrasound system.
5. Clinical image acquired by the HS60 V2.03 ultrasound system.
6. S-Detect™ for Breast and S-Detect™ for Thyroid are not available in Canada.
7. In the United States, certain classification items in S-Detect™ such as Margin, Posterior Features, and Echo Pattern are not automatically provided. Also the recommendations about whether results are benign or malignant in S-Detect™ are not applicable in the United States.
8. Strain value for ElastoScan+™ is not applicable in Canada and the United States.
9. A purchase of Mobile Export option is required to use HelloMom™.

SAMSUNG MEDISON CO., LTD.

© 2021 Samsung Medison All Rights Reserved.
Samsung Medison reserves the right to modify the design, packaging, specifications, and features shown herein, without prior notice or obligation.

Focus on your needs

Ultrasound system HS60 Powered by CrystalLive™



CT-HS60 V2.04-OB-IMC-211206-EN

Scan code or visit
www.samsunghealthcare.com
to learn more



Life-long healthcare solutions for women

Samsung is committed to supporting life-long healthcare for women by diligently pursuing innovative healthcare solutions in Obstetrics and Gynecology ultrasound.

Samsung's HS60 ultrasound system adopts the integrated solution, combining high-quality imaging with advanced features to support healthcare professionals in making fast and accurate decisions.



Powered by CrystalLive™


CrystalLive™ is Samsung's ultrasound imaging engine with enhanced image processing capabilities including 2D image processing, 3D rendering and color signal processing. It provides realistic and coherent images in various modes. In addition, it also offers an efficient diagnostic environment by handling large amount of signal data quickly.

Samsung
Ultrasound System **HS60**




S-Vue Transducer™

S-Vue Transducer™ provides more efficient piezoelectric properties, resulting in wider bandwidths that enable better penetration and higher quality resolution.



S-Vue Transducer™

CA1-7AD, CA2-9AD, CA3-10A,
CV1-8AD, PA1-5A



* The image is for illustrational purposes only and might differ from the actual performance of the device.

Focusing on diagnostic solutions for women's health

The HS60 provides a versatile range of obstetric and gynecological capabilities for efficient and effective care. With its efficient diagnostic solutions, the HS60 supports your knowledge and experience to help you make clear, confident decisions.



Family Planning

Measure the size of follicles based on 2D

2D Follicle™¹ is a function to measure the size of follicles based on 2D image and to provide information about the status during gynecologic examination.

Assess the risk of infertility


5D Follicle™¹ measures the size of multiple follicles and provides information on their status during OB exams.



Healthy Pregnancy Biometry

Measure fetal biometry parameter in one click

A semi-automatic technology for biometric measurement, **BiometryAssist™**¹, enables users to measure the growth of the fetus more quickly and with greater accuracy while maintaining exam consistency.



Fetal brain³




Healthy Pregnancy Diagnosis

Measure NT using automatic detection of mid-sagittal plane

5D NT™¹ provides the midsagittal plane view automatically by rotating and magnifying the images when measuring the nuchal translucency (NT) of the fetus in early weeks.

Examine fetal heart including blood flow dynamics

5D Heart Color™¹ provides 9 standard planes of the heart by using the fetal STIC data as well as important information about fetal heart development in an easy and accurate way in accordance with the AIUM guideline. In addition, it offers dedicated Preset, Predictive Cursor, Diagnostic Alert, and heart Diastole/Systole timepoints.




Nuchal translucency²



Healthy Pregnancy Visualization

Visualize internal and external structures by volume rendering

CrystalVue™¹ is an advanced volume rendering technology that enhances visualization of both internal and external structures in a single rendered image using a combination of intensity, gradient and position.



Fetus³ in 3D

Express 3D anatomy in detailed and realistic view

RealisticVue™¹ displays high resolution 3D anatomy with detailed expression and realistic depth perception. User selectable light source direction creates intricately graduated shadows for better defined anatomical structures.



Fetal face³ with RealisticVue™




Healthy Birth

Support in deciding delivery method

LaborAssist™¹ is a function that provides information of the progress of delivery by the automatic measurement of AoP (Angle of Progression) and the direction of the fetal head. This not only helps in effective communication between the healthcare professionals and mothers, but also assists in deciding a delivery method for the healthcare professionals.

* AoP complies with the metrics specified in the ISUOG Guideline.



AoP measurement⁴ with LaborAssist™



Gynecology & Breast Health

Analyze selected breast lesions and report breast assessment

S-Detect™ for Breast^{1,6,7} analyzes selected lesions in the breast ultrasound study and shows the analysis data, applies BI-RADS ATLAS* (Breast Imaging-Reporting and Data System, Atlas) to provide standardized reporting; and helps diagnosis with the streamlined workflow.

* It is a registered trademark of ACR and all rights reserved by ACR.




Breast assessment³ with S-Detect™ for Breast

Highly detailed images through innovation


Samsung's innovative imaging technologies and S-Vue Transducer™ provide highly detailed images to increase diagnostic confidence.

Reduce noise to improve 2D image quality

The noise reduction filter improves edge enhancement and creates sharper 2D images for optimal diagnostic performance. In addition, **ClearVision** provides application-specific optimization and advanced temporal resolution in live scan mode.




Fetal brain²



Fetal brain² with **ClearVision**

Visualize slow flow in microvascular structures

MV-Flow™¹ visualizes microcirculatory and slow blood flow to display the intensity in color. It is suitable for observation of microcirculatory and volume of slow blood flow.




Fetal Brain⁵ with MV-Flow™




Umbilical Cord⁵ (MV-Flow™ with LumiFlow™)

Show blood flow in vessels in a 3D like display


LumiFlow™¹ is a function that visualizes blood flow in three dimensional-like to help understand the structure of blood flow and small vessels intuitively.




Fetal brain³ with S-Flow™



Umbilical cord³ with S-Flow™




Fetal face³ with Mirror View




NT³ measured with 2D NT




Fetal brain³ with MSV



Fetal heart² with PW



Fetal heart² with color



4 chamber view²


Enhanced Efficiency

The HS60 has been designed to enhance efficiency through reducing keystrokes, enabling you to streamline your workflow by combining multiple actions into one. Its user-oriented design also enables you to focus on your patient, reducing the complexity and stress of operating the system.



Build predefined protocols for streamlined process

EzExam+™ enables you to build or use a predefined protocol, and assign protocols for examinations that are regularly performed in the hospital in order to reduce the number of steps that you have to go through. For fetus diagnosis, in particular, you can arrange the examination order according to the fetus position using the touchscreen, and automatically apply the BodyMarker, Annotation, Measurement, etc.




Compare previous and current exam in a side-by-side display

EzCompare™ allows easy access to previously taken exams to evaluate corresponding views in a side-by-side display. For greater efficiency, EzCompare™ automatically matches the image settings, annotations, and bodymarkers from the prior study.




Magnify the region of interest in a picture-in-picture window

When placing a caliper, **Measure Navigation** automatically magnifies the area of interest using a picture-in-picture window to allow more precise placement of the calipers. This is especially useful when measuring small structures or when accuracy is critical.




Select transducer and preset combinations in one click

With one touch, the user can select the most common transducer and preset combinations. **QuickPreset** increases efficiency to make a full day of scanning simple and easy.



View images in wider screen

WideScreen provides approximately 27% more lateral viewing information compared to normal screen, allowing ultrasonic examination with wider view at a glance.




Simple transfer of fetal ultrasound images and clips

HelloMom™ 1.9 is a simple and secure image sharing solution that generates a QR code for the selected fetal images to be transferred. HelloMom™ allows pregnant women and their family to download fetal ultrasound images simply by scanning the QR code with their smartphones, reducing the hassle of installing a separate application.



1 Gel warmer¹

Two-level adjustable gel warmer maintains ultrasound gel at a comfortable temperature.



3 Use the system when AC power is temporarily unavailable

BatteryAssist™ 1 provides the system with battery power. It enables users to perform scans and transport the ultrasound system to other locations in environments where AC power may not be available temporarily.




2 Solid State Drive (SSD)

The HS60 uses advanced solid state drives. These stable and dependable drives allow faster bootup, better frame rates, and fast processing speeds.

4 Clever use of space

With its reduced weight and compact size, the HS60 takes up minimal space and can move freely. In addition, its streamlined rear profile allows you to park the HS60 in small spaces.



Comprehensive selection of transducers

Curved array transducers



CA1-7AD
abdomen, obstetrics,
gynecology,
musculoskeletal, pediatric,
vascular, urology



CA2-9AD
abdomen, obstetrics,
gynecology,
musculoskeletal, pediatric,
vascular, urology



CA3-10A
abdomen, obstetrics,
gynecology,
musculoskeletal, pediatric,
vascular, urology



CA4-10M
abdomen, obstetrics,
gynecology,
musculoskeletal, pediatric,
vascular, urology

Linear array transducers



LA3-14AD
small parts, vascular,
musculoskeletal,
abdomen, obstetrics,
gynecology, pediatric



LA3-16A
small parts, vascular,
musculoskeletal,
abdomen, obstetrics,
gynecology, pediatric



LA2-9A
small parts, vascular,
musculoskeletal,
abdomen, obstetrics,
gynecology, pediatric



LA4-18BD
small parts, vascular,
musculoskeletal,
abdomen, obstetrics,
gynecology, pediatric



LA3-16AI
small parts, vascular,
musculoskeletal,
abdomen, obstetrics,
gynecology, pediatric

Volume transducers



CV1-8AD
abdomen, obstetrics,
gynecology,
musculoskeletal, pediatric,
vascular, urology



EV2-10A
obstetrics, gynecology,
urology

Endo-cavity transducers



EA2-11AV
obstetrics, gynecology,
urology



EA2-11AR
obstetrics, gynecology,
urology



EA2-11B
obstetrics, gynecology,
urology

Phased array transducers



PA1-5A
abdomen, cardiac,
vascular, pediatric



PA3-8B
abdomen, cardiac,
vascular, pediatric



PA4-12B
abdomen, cardiac,
vascular, pediatric

CW transducers



CW6.0
cardiac, vascular



DP2B
cardiac, vascular

TEE transducers



MMPT3-7
cardiac

* Some of the transducers may not be available in some countries.